

Surgical Instruments Catalogue
Chirurgische Instrumente Katalog
Instruments chirurgicaux catalogue

DERMATOLOGY INSTRUMENTS



Manufacturer & Exporters of Surgical Instruments



**SHEFFIELD
SURGICAL
INSTRUMENTS**

Quality

Instruments

Since 1992



Quality Instruments Since 1992
Manufacturer & Exporters of Surgical Instruments

Sheffield Surgical Building, Wazirabad Road, Ugoki, Sialkot-51310 Pakistan.

Tell # +92-333-1400682

Email: sheffieldsurgicalinstruments@gmail.com


Email: info@sheffieldsurgicalinst.com

Web: www.sheffieldsurgicalinst.com



Dermatology

SCISSORS • SUPERCUT



SuperCut Scissors are manufactured with a razor edge on one blade and micro serrations on the other blade. The razor edge blade is designed to provide a smooth, precise and less traumatic cutting of tissue and the micro serrated blade is designed to reduce slippage of tissue within the blades of the scissor. The SuperCut is only recommended for tissue not for suture or gauze. SuperCut Scissors are designated by Black Ring Handles.

IRIS SCISSORS • SC

- SSI-08-34011SC IRIS Scissor, straight, 4½" SuperCut
- SSI-08-34111SC IRIS Scissor, curved, 4½" SuperCut



METZENBAUM SCISSORS • SC

- SSI-08-28015SC METZENBAUM Scissor, straight, 5¾" SuperCut
- SSI-08-28115SC METZENBAUM Scissor, curved, 5¾" SuperCut



STRABISMUS SCISSORS • SC

- SSI-08-37011SC STRABISMUS/Baby Metzenbaum Scissor, straight, 4¼" SuperCut
- SSI-08-37111SC STRABISMUS/Baby Metzenbaum Scissor, curved, 4¼" SuperCut




STEVENS TENOTOMY SCISSORS • SC

- SSI-08-36011SCS STEVENS Tenotomy Scissor, straight, sharp, 4¼" SuperCut
- SSI-08-36111SCS STEVENS Tenotomy Scissor, curved, sharp, 4¼" SuperCut





SCISSORS • TUNGSTEN CARBIDE



Tungsten Carbide (TC) Scissors are manufactured with Tungsten Carbide Inserts in the blades of the Scissor. Tungsten Carbide is much harder and longer lasting than Stainless Steel and is designed to provide a longer lasting sharp edge on the Scissor blade for many years. Tungsten Carbide Scissors are designated by Gold Ring Handles.

IRIS SCISSORS • TC

- SSI-08-34411 IRIS Scissor, straight, 4½" TC
- SSI-08-34511 IRIS Scissor, curved, 4½" TC



OR SCISSORS • TC

- SSI-08-12413 OR Scissor, straight, sharp/sharp, 5" TC
- SSI-08-12414 OR Scissor, straight, sharp/sharp, 5¾" TC



METZENBAUM SCISSORS • TC

- SSI-08-28415 METZENBAUM Scissor, straight, blunt, 5¾" TC
- SSI-08-28515 METZENBAUM Scissor, curved, blunt, 5¾" TC



STRABISMUS SCISSORS • TC

- SSI-08-37411 STRABISMUS/Baby Metzenbaum Scissor, straight, 4¼" TC
- SSI-08-37511 STRABISMUS/Baby Metzenbaum Scissor, curved, 4¼" TC



KILNER SCISSOR • TC

- SSI-08-31412 KILNER Undermining Scissor, curved, non-serrated, flat tips, 4¾"



Dermatology

SCISSORS

IRIS SCISSORS

- SSI-08-33810 IRIS Scissor, straight, 4"
- SSI-08-33910 IRIS Scissor, curved, 4"
- SSI-08-34011 IRIS Scissor, straight, 4½"
- SSI-08-34111 IRIS Scissor, curved, 4½"



IRIS RIBBON SCISSORS

- SSI-08-32209 IRIS Scissor, round, straight, 17mm blade, ribbon type, 3¾"
- SSI-08-32309 IRIS Scissor, round, curved, 17mm blade, ribbon type, 3¾"
- SSI-08-32610 IRIS Scissor, round, straight, 22mm blade, ribbon type, 4⅞"
- SSI-08-32710 IRIS Scissor, round, curved, 22mm blade, ribbon type, 4⅞"



OR SCISSORS

- SSI-08-12011 OR Scissor, straight, sharp/sharp, 4½"
- SSI-08-12014 OR Scissor, straight, sharp/sharp, 5"



METZENBAUM SCISSORS

- SSI-08-28015 METZENBAUM Scissor, straight, blunt, 5¾"
- SSI-08-28115 METZENBAUM Scissor, curved, blunt, 5¾"



STRABISMUS SCISSORS

- SSI-08-37011 STRABISMUS/Baby Metzenbaum Scissor, straight, 4¼"
- SSI-08-37111 STRABISMUS/Baby Metzenbaum Scissor, curved, 4¼"



KILNER SCISSOR

- SSI-08-31012 KILNER Undermining Scissor, curved, non-serrated, flat tips, 4¾"





Dermatology

LITTLER DISSECTING SCISSOR

SSI-08-20512 LITTLER Dissecting Scissor, curved tips with eye for suture, 4¾"



STEVENS TENOTOMY SCISSORS

SSI-08-36011 STEVENS Tenotomy Scissor, straight, sharp, 4¼"

SSI-08-36111 STEVENS Tenotomy Scissor, curved, sharp, 4¼"



GRADLE SCISSOR

SSI-08-34653 GRADLE Eye Suture Scissor, ribbon handle, slightly curved, sharp, 3¾"



SPENCER STITCH SCISSORS

SSI-08-80209 SPENCER Stitch Scissor, 3½"

SSI-08-80210 SPENCER Stitch Scissor, 4¼"



STITCH SCISSOR

SSI-08-84664 Stitch Scissor, angled, 4½"



WILMER CONVERSE SCISSOR

SSI-08-35010 WILMER CONVERSE Scissor, angled on flat, 4"



NOYES IRIS SCISSORS


SSI-08-11012 NOYES IRIS Scissor, spring handle, straight, sharp, 4¾"

SSI-08-11112 NOYES IRIS Scissor, spring handle, curved, sharp, 4¾"



Dermatology

NEEDLE HOLDERS • TUNGSTEN CARBIDE



Tungsten Carbide (TC) Needle Holders are manufactured with Tungsten Carbide Inserts in the jaws of the Needle Holder. Tungsten Carbide is much harder and longer lasting than Stainless Steel and provides sharper teeth which are designed to provide a more precise and secure grip on the needles. Tungsten Carbide Needle Holders are designated by Gold Ring Handles.

WEBSTER NEEDLE HOLDERS • TC

- SSI-24-104429 WEBSTER Needle Holder, cross serrated jaw, 5" TC
- SSI-24-104490 WEBSTER Needle Holder, smooth, 5" TC



MAYO-HEGAR NEEDLE HOLDERS • TC

- SSI-24-19414 MAYO-HEGAR Needle Holder, 5" TC
- SSI-24-19416 MAYO-HEGAR Needle Holder, 6" TC




DERF NEEDLE HOLDER • TC

- SSI-24-11412 DERF Needle Holder, 4¾" TC



CRILE-WOOD NEEDLE HOLDER • TC

- SSI-24-16415 CRILE-WOOD Needle Holder, 6" TC



OLSEN-HEGAR NEEDLE HOLDERS • TC

- SSI-24-41412 OLSEN-HEGAR Needle Holder with Suture Scissor, serrated, 4¾" TC
- SSI-24-41414 OLSEN-HEGAR Needle Holder with Suture Scissor, serrated, 5½" TC





NEEDLE HOLDERS

WEBSTER NEEDLE HOLDER


SSI-24-104480 WEBSTER Needle Holder, smooth, 5"



MAYO-HEGAR NEEDLE HOLDERS

SSI-24-18014 MAYO-HEGAR Needle Holder, 5"

SSI-24-18016 MAYO-HEGAR Needle Holder, 6"



DERF NEEDLE HOLDER

SSI-24-11012 DERF Needle Holder, 4¾"



CRILE-WOOD NEEDLE HOLDER

SSI-24-16015 CRILE-WOOD Needle Holder, 6"



OLSEN-HEGAR NEEDLE HOLDERS


SSI-24-41012 OLSEN-HEGAR Needle Holder with Suture Scissor, straight, serrated, 4¾"

SSI-24-41014 OLSEN-HEGAR Needle Holder with Suture Scissor, straight, serrated, 5½"



Dermatology


TISSUE FORCEPS • TUNGSTEN CARBIDE



Tungsten Carbide (TC) Tissue Forceps are manufactured with Tungsten Carbide Inserts in the tips of the Forcep. The TC platform is in front of the 1x2 teeth on the distal tip and is used to provide a secure grip of the needle, suture or tissue. Tungsten Carbide is much harder and longer lasting than Stainless Steel. Tungsten Carbide Tissue Forceps are designated by Gold on the proximal end of the Forcep.

ADSON FORCEPS


- SSI-10-17412 ADSON Forcep, serrated, 4¾" TC
- SSI-10-17512 ADSON Forcep, 1x2 teeth, with tying platform, 4¾" TC



TISSUE FORCEPS

ADSON FORCEPS

- SSI-10-17012 ADSON Dressing Forcep, delicate, serrated, 4¾"
- SSI-10-18012 ADSON Tissue Forcep, delicate, 1x2 teeth, 4¾"
- SSI-10-180941 ADSON TUBIANA Tissue & Suture Forcep, tying platform, 1x2 teeth, 4¾"
- SSI-10-182609 ADSON Forcep, cross serrations, 1x2 teeth, 4¾"
- SSI-10-18613 ADSON Tissue Forcep, fenestrated handle, delicate, 1x2 teeth, 4¾"
- SSI-10-18614 ADSON Tissue Forcep, fenestrated handle, delicate, serrated, 4¾"
- SSI-10-18812 BROWN-ADSON Tissue Forcep, 7x7 side-grasp teeth, 4¾"



BISHOP-HARMON FORCEPS

- SSI-42-94005 BISHOP HARMON Iris Forcep, straight, serrated, 0.5mm tip, 3½"
- SSI-42-94008 BISHOP HARMON Iris Forcep, straight, serrated, 0.8mm tip, 3½"
- SSI-42-94206 BISHOP HARMON Iris Forcep, straight, 1x2 teeth, 0.3mm tip, 3½"
- SSI-42-94205 BISHOP HARMON Iris Forcep, straight, 1x2 teeth, 0.5mm tip, 3½"
- SSI-42-94208 BISHOP HARMON Iris Forcep, straight, 1x2 teeth, 0.8mm tip, 3½"



CILIA FORCEP

- SSI-43-42710 Cilia and Suture Forcep, precision fitted diagonal jaws 3mm wide, 3¾"



DESMARRES CHALAZION FORCEPS

- SSI-43-25409 DESMARRES Chalazion Forcep, small, 3½"
- SSI-43-25609 DESMARRES Chalazion Forcep, medium, 3½"
- SSI-43-25809 DESMARRES Chalazion Forcep, large, 3½"





HEMOSTATIC FORCEPS

HARTMAN HEMOSTATIC FORCEPS

- SSI-12-21010 HARTMAN Hemostatic Forcep, straight, 4"
SSI-12-21110 HARTMAN Hemostatic Forcep, curved, 4"



HALSTED HEMOSTATIC FORCEPS

- SSI-12-22012 HALSTED Hemostatic Forcep, straight, 5"
SSI-12-22112 HALSTED Hemostatic Forcep, curved, 5"



BACKHAUS FORCEPS

- SSI-14-11109 BACKHAUS Towel Clamp, 3½"
SSI-14-11113 BACKHAUS Towel Clamp, 5¼"



BLADE HANDLES

BLADE HANDLE

- SSI-06-10301 Blade Handle No. 3, graduated



SIEGEL BLADE HANDLE

- SSI-06-10303 SIEGEL Round Blade Handle



BEAVER BLADE HANDLES

- SSI-06-10910 BEAVER Scalpel Blade Handle, 4"
SSI-06-10913 BEAVER Scalpel Blade Handle, 5"



Dermatology

CURETTES



FOX CURETTES

- SSI-75-12801 FOX Curette (Ø 1mm), oval, 5½"
- SSI-75-12802 FOX Curette (Ø 2mm), oval, 5½"
- SSI-75-12803 FOX Curette (Ø 3mm), oval, 5½"
- SSI-75-12804 FOX Curette (Ø 4mm), oval, 5½"
- SSI-75-12805 FOX Curette (Ø 5mm), oval, 5½"
- SSI-75-12806 FOX Curette (Ø 6mm), oval, 5½"

- SSI-75-12901 FOX Curette (Ø 1mm), round, 5½"
- SSI-75-12902 FOX Curette (Ø 2mm), round, 5½"
- SSI-75-12903 FOX Curette (Ø 3mm), round, 5½"
- SSI-75-12904 FOX Curette (Ø 4mm), round, 5½"
- SSI-75-12905 FOX Curette (Ø 5mm), round, 5½"
- SSI-75-12906 FOX Curette (Ø 6mm), round, 5½"

EXTRACTORS

SCHAMBERG EXTRACTORS

- SSI-74-11101 SCHAMBERG Comedone Extractor, slightly curved double-ends, square/crimped small loop, 3¾"
- SSI-74-11102 SCHAMBERG Comedone Extractor, slightly curved double-ends, extra fine loops, square/crimped small loop, 3¾"
- SSI-74-11103 SCHAMBERG Comedone Extractor, slightly curved double crimped ends, 3¾"
- SSI-74-11110 SCHAMBERG Comedone Extractor, slightly curved double square ends, 3¾"
- SSI-74-11111 SCHAMBERG Comedone Extractor, double square ends, 3¾"



WALTON COMEDONE EXTRACTOR

- SSI-74-11814 WALTON Comedone Extractor, double-ended, sharp/round, 5½"



UNNA COMEDONE EXTRACTORS

- SSI-74-11314 UNNA Comedone Extractor, spoon ends, 5½"
- SSI-74-11514 UNNA Comedone Extractor, round ends, 5½"



SAALFELD COMEDONE EXTRACTOR

- SSI-74-12110 SAALFELD Comedone Extractor, double-ended, sharp/round, 4"



ELEVATORS, RETRACTORS & SKIN HOOKS

LOCKE ELEVATOR

- SSI-74-74011 Locke Elevator, small blade, 3x15mm, 4½"



DESMARRES RETRACTORS

- SSI-18-23010 DESMARRES Eyelid Retractor, 10mm blade, 6¼"
- SSI-18-23012 DESMARRES Eyelid Retractor, 12mm blade, 6¼"



CROSS ACTION RETRACTORS

- SSI-18-63410 Cross Action Retractor, sharp, 4"
- SSI-18-63510 Cross Action Retractor, blunt, 4"





Dermatology

COTTLE SKIN HOOK

SSI-46-55700 COTTLE Type Skin Hook, angled, 5½"



FRAZIER SKIN HOOK

SSI-18-21413 FRAZIER Dura Hook, sharp, 1 prong, 5½"



SSI- SKIN

SSI-18-20001 Skin Hook, single prong, sharp, 6¼"

SSI-18-20002 Skin Hook, double prong, 5mm, sharp, 6¼"



GUTHRIE SKIN HOOKS

SSI-18-21201 GUTHRIE Double Hook, small, 2.5mm, 4¾"

SSI-18-21202 GUTHRIE Double Hook, large, 3mm, 4¾"



JOSEPH SKIN HOOKS


SSI-18-21601 JOSEPH Skin Hook, 1 prong, 6¼"

SSI-18-21602 JOSEPH Skin Hook, 2 prongs, 2mm, 6¼"

SSI-18-21605 JOSEPH Skin Hook, 2 prongs, 5mm, 6¼"

SSI-18-21607 JOSEPH Skin Hook, 2 prongs, 7mm, 6¼"

SSI-18-21610 JOSEPH Skin Hook, 2 prongs, 10mm, 6¼"



NAIL NIPPERS

NAIL NIPPER

SSI-74-33014 Nail Nipper, straight, double spring, with lock 5½"



NAIL NIPPERS • BARREL SPRINGS

SSI-74-33615 Nail Nipper, double action, sharp concave jaw, barrel spring, 6"

SSI-74-33616 Nail Nipper, double action, sharp straight jaw, barrel spring, 6"



ENGLISH ANVIL NAIL SPLITTER

SSI-74-33413 English Anvil Pattern Nail Splitter, double spring, 5"





Sheffield Surgical Building, Wazirabad Road, Ugoki, Sialkot-51310 Pakistan.

Tell # +92-333-1400682

Email: sheffieldsurgicalinstruments@gmail.com

Email: info@sheffieldsurgicalinst.com

Web: www.sheffieldsurgicalinst.com

