

Surgical Instruments Catalogue  
Chirurgische Instrumente Katalog  
Instruments chirurgicaux catalogue

## OPHTHALMIC INSTRUMENTS



Manufacturer & Exporters of Surgical Instruments



**SHEFFIELD  
SURGICAL  
INSTRUMENTS**

Quality

Instruments

Since 1992



**Quality**

**Instruments**

**Since 1992**

**Manufacturer & Exporters of Surgical Instruments**

**Sheffield Surgical Building, Wazirabad Road, Ugoki, Sialkot-51310 Pakistan.**

**Tell # +92-333-1400682**

**Email: sheffieldsurgicalinstruments@gmail.com**

**Email: info@sheffieldsurgicalinst.com**

**Web: www.sheffieldsurgicalinst.com**





**Gill**  
Iris Forceps, curved  
criss-cross serrated tips  
1x2 teeth  
FC-2001



**Botvin**  
Iris Forceps  
for peripheral iridectomy  
serrated tips  
1x2 teeth  
FC-2002



**Gills-Colibri**  
Forcep  
very fine pointed tips  
FC-2003



**Colibri**  
Forcep, very delicate  
1x2 teeth, 0.12mm  
FC-2004



**Harms-Colibri**  
Forcep, very delicate  
1x2 teeth, 0.12mm  
with tying platform  
FC-2005



**Dyson-Colibri**  
Forcep, very delicate  
1x2 teeth, 0.12mm  
extremely thin shanks  
short tips with tying platform  
FC-2006



# Forceps

## Colibri

Forceps, standard  
1x2 teeth, 0.4mm  
with tying platform  
FC-2007



## Polack

Corneal suturing forceps  
1mm spread, very delicate  
1x2 teeth, 0.12mm  
FC-2008



## Hofmann-polack

Corneal suturing forceps  
with guide pin, 1mm spread  
1x1pierce teeth, 0.2mm  
FC-2009




## Pierse-colibri

Corneal forceps  
very delicate  
1x1pierce teeth, 0.1mm  
FC-2010



## Pierse

Corneal forceps  
very delicate  
1x1pierce teeth, 0.1mm  
straight  
curved  
FC-2011



**Castroviejo-Colibri**  
Forceps very delicate  
1x2 teeth, 0.12mm  
FC-2012





**Girards**  
 Corneoscleral Forceps  
 1x2 teeth, 0.12mm  
 with tying platform  
 right handed  
 left handed  
 FC-2013



**Maumenee-Colibri**  
 Corneoscleral Forceps  
 very delicate  
 1x2 teeth, 0.12mm  
 FC-2014



**Tennant-Colibri**  
 Corneoscleral Forceps  
 very delicate  
 1x2 teeth, 0.12mm  
 FC-2015



**Bores**  
 Corneoscleral Fixation Forceps  
 U-shaped, 3mm spread  
 1x2 teeth at each tip,  
 0.12mm  
 FC-2016



**Rubman Bores**  
 Corneoscleral Fixation Forceps  
 U-shaped, 3mm spread  
 1x2 teeth at each tip,  
 0.12mm  
 FC-2017



**Arrowsmith**  
 Corneoscleral Fixation Forceps  
 with lock  
 U-shaped, 8mm spread  
 1x2 teeth at each tip,  
 0.12mm  
 FC-2018S  
 FC-2018A



# Forceps



## Dressing

Forceps, serrated  
delicate  
straight  
curved

FC-2024S  
FC-2014A



## Tissue

Tissue Forceps  
1x2 teeth  
delicate  
standard

FC-2025S  
FC-2025A



## Tissue

Forceps, 1x2 teeth  
delicate  
straight  
curved

FC-2026S  
FC-2026A



## Lester

Fixation Forceps  
straight  
1x2teeth  
2x3teeth

FC-2027S  
FC-2027A



## Barken

Fixation Forceps  
1x2 teeth

FC-2028





**Katena**  
 Corneal Fixation Forceps  
 with lock  
 U-shaped, 13mm spread  
 1x2 teeth at each tip,  
 0.12mm  
 straight  
 angled  
**FC-2019S**  
**FC-2019A**



**Bishop-Harmon**  
 Tissue Forceps  
 1x2 teeth  
 delicate  
 standard  
**FC-2020S**  
**FC-2020T**



**Bishop-Harmon**  
 Dressing Forceps  
 criss-cross serrated tips  
**FC-2021**



**Sauer**  
 Fixation Forceps  
 1x2 angled teeth  
**FC-2022**



**Manhattan Eye & Ear**  
 Fixation Forceps  
 1x2 angled teeth  
 with tying platform  
**FC-2023**



# Forceps

**Farkas-Bracken**  
Fixation Forceps  
1x2teeth, delicate  
FC-2029



**Kuhnt**  
Fixation Forceps  
1x2 curved teeth  
FC-2030



**Elschnig**  
Fixation Forceps  
1x2 angled teeth  
FC-2031



**Nugent**  
Superior Rectus Forceps  
angled, 1x1 teeth  
FC-2032



**Troutman**  
Superior Rectus Forceps  
angled, 1x2 teeth  
FC-2033



**Tennant-Troutman**  
Superior Rectus Forceps  
curved shanks, 1x2 teeth  
FC-2034



**Graefe**  
Fixation Forceps  
4.5mm wide jaws  
without lock  
with lock  
FC-2035A  
FC-2035B






# Forceps



## Tennant

Tying Forceps  
with guide pin  
extra delicate  
smooth jaws with  
6mm long platform  
for 9-0 to 11-0 sutures  
straight  
curved



**FC-2036S**  
**FC-2036C**

## Harms

Tying Forceps  
smooth jaws  
for 7-0 to 9-0 sutures  
straight  
curved



**FC-2037S**  
**FC-2037C**

## Fine

Tying Forceps  
smooth curved jaws  
for 7-0 to 9-0 sutures  
**FC-2038**



## Barraquer

Cilia Forceps  
lightly angled  
smooth curved jaws  
**FC-2039**



## Katena

Cilia Forceps  
slanted jaws  
gold plated  
**FC-2040**



## Gradle

Cilia Forceps  
smooth, squared off  
positive gripping jaws  
**FC-2041**



# Forceps

**Lambert**  
Chalazion Forceps  
round, 8mm diameter  
FC-2042



**Hunt**  
Chalazion Forceps  
round, 8mm diameter  
FC-2043



**Baird**  
Chalazion Forceps  
oval, small  
FC-2044



**Francis**  
Chalazion Forceps  
oval, medium  
FC-2045



**Desmarres**  
Chalazion Forceps  
oval, large  
FC-2046



**Wies**  
Chalazion Forceps  
grooved upper jaw  
FC-2047



**Lordan**  
Chalazion Forceps  
triangular jaws  
FC-2055



# Forceps



**Jewelers  
Forceps, delicate  
FC-2048**



**Jewelers  
Forceps, delicate  
FC-2049**



**Jewelers  
Forceps, extra delicate  
FC-2050**



**Osher  
Iris Forceps  
very fine textured  
jaws  
FC-2051**



**Jewelers  
Forceps, extra delicate angled  
FC-2052**



**Jewelers  
Forceps, curved  
FC-2053**



**Suture Removal Forceps  
Very fine pointed tips  
with carbide impregnated jaws  
FC-2054**



# Forceps

**Serrefine**  
Small, straight  
FC-2055



**Hartman**  
Mosquito Forceps  
extra delicate  
3 3/4"(95mm) long  
straight  
curved  
FC-2056



**Halsted**  
Mosquito Forceps  
delicate  
5 3/4"(145mm) long  
straight  
curved  
FC-2057



**Backhaus**  
Towel Clamp  
3 1/2"(90mm)long  
5 1/4"(135mm) long  
FC-2058



**Baby Jones**  
Towel Clamp  
cross-action  
2 1/4"(55mm) long  
FC-2059



# Forceps

**Shepard II**  
 Lens Holding Forceps  
 with nudge notch  
 for A/C IoLs  
 with lock  
 without lock  
**FC-2060**



**Jameson**  
 Muscle Forceps  
 with slide lock  
 4teeth, child size  
 right  
 left  
**FC-2061R**  
**FC-2061L**



**Jameson**  
 Muscle Forceps  
 with slide lock  
 6teeth, adult size  
 right  
 left  
**FC-2062R**  
**FC-2062L**



**Berke**  
 Ptosis Forceps  
 with slide lock  
 longitudinally grooved jaws  
 20mm long jaws  
 27mm long jaws  
**FC-2063A**  
**FC-2063B**



**Erhardt**  
 Lid Forceps  
**FC-2064**



# Scissors

**Castroviejo**  
Keratoplasty Scissors  
angled to side  
blunt tips  
small blades  
SC-2065



**Castroviejo**  
Keratoplasty Scissors  
angled to side  
blunt tips  
medium blades  
SC-2066



**McPherson-Westcott**  
Conjunctival Scissors  
curved, blunt tips  
small blades  
SC-2067



**Westcott**  
Tenotomy Scissors  
curved, blunt tips  
medium blades  
SC-2068



**Westcott**  
Tenotomy Scissors  
curved, blunt tips  
standards blades  
right  
left  
SC-2069R  
SC-2069L



# Markers & Calipers




## Anis

Radial Marker  
40mm center opening  
12.5mm O.D.

- |          |           |
|----------|-----------|
| 3 blades | 10 blades |
| 4 blades | 12 blades |
| 6 blades | 16 blades |
| 8 blades |           |

MC-2070A

MC-2070B



## Ains

Suture Placement Marker  
for corneal grafts, 8 radial blades  
with 8mm center opening

MC-2071



## Deitz

Incision Depth Gauge  
for corneal incisions

0.5mm

0.55mm

MC-2072



## Thornton

Limbal Incision Ruler  
with marking lines  
lines spaced 1mm apart\*

MC-2073



## Lindstrom

Small Incision Marker  
for cataract incision widths of  
3, 4, and 5mm at 1, 2, or 3mm  
away from the limbus\*


MC-2074



## Castroviejo

Caliper, 20mm spread  
straight  
curved

MC-2075



## Osher

Incision Caliper  
for internal measurement of  
cataract incisions from 1 to 10mm

MC-2076



# Scissors



**Castroviejo**  
**Corneoscleral Scissors**  
 small blades  
 right  
 left  
 MC-2077R  
 MC-2077L



**Castroviejo**  
**Corneoscleral Scissors**  
 Medium blades  
 Right  
 Left  
 MC-2078R  
 MC-2078L



**Troutman-Castroviejo**  
**Corneoscleral Scissors**  
 with stop, small blades  
 Right  
 Left  
 MC-2079R  
 MC-2079L





# Scissors

**Vannas**  
Capsulotomy Scissors  
sharp tips  
angled to side  
SC-2080



**Gills-Vannas**  
Capsulotomy Scissors  
extra thin long blades  
sharp tips  
medium blades  
SC-2081



**Gills-Vannas**  
Capsulotomy Scissors  
extra thin long blades  
sharp tips  
SC-2082



**Gills-Vannas**  
Capsulotomy Scissors  
extra thin long blades  
sharp tips  
SC-2083



**McPherson-Vannas**  
Iris Scissors  
sharp tips  
curved  
SC-2084



**Katena-Vannas**  
Scissors, extra thin  
7mm long blades  
sharp points  
straight  
SC-2085



# Forceps

## Shaepard

Lens holding Forceps  
with serrated lower jaw  
gently curved shanks  
without lock  
with lock  
FC-2086



## Sinsky

Lens holding Forceps  
delicate T-shaped lower  
jaw for 3-points fixation  
FC-2087



## Kraff

Lens holding Forceps  
delicate curved jaws  
to grasp edge of lens optic  
FC-2088



## Anis

Lens Inserting Forceps  
very delicate curved tips  
to grasp IOL haptics.  
Angled shanks for  
comfortable hand position  
FC-2089



## Clayman

Lens holding Forceps  
delicate angled jaws  
gently curved shanks  
without lock  
with lock  
FC-2090



## Phaco-Folder II

For folding and inserting  
soft IOLs through a  
small incision  
FC-2091



# Scissors

## Eye Scissors

Standard  
pointed tips  
straight  
curved  
SC-2092S  
SC-2092C




## Eye Scissors

Standard, large rings  
pointed tips  
straight  
curved  
SC-2093S  
SC-2093C



## Stitch Scissors

Ribbon type  
needle points  
light curved  
SC-2094



# Spatulas

**Barraquer**  
 Corneoscleral brush  
 Sable with plastic handle  
 SP-2095



**McPherson**  
 Iris Spatula  
 1mm wide  
 malleable stainless steel  
 SP-2096



**Culler**  
 Iris spatula  
 Sterling silver  
 Blade  
 1mm wide  
 2mm wide  
 SP-2097A  
 SP-2097B



**Wheeler**  
 Double-Ended Spatula  
 Lightly curved  
 SP-2098



**Barraquer**  
 Iris Spatula  
 Round, smooth  
 0.5mm dia.  
 0.25 mm dia.  
 SP-2099



**Gimbel**  
 Nucleus Spatula  
 Angled, round, smooth  
 0.25mm dia.  
 SP-2100



**Knolle**  
 Lens Nucleus Spatula  
 Malleable stainless steel  
 Flat, 0.5mm wide  
 SP-2101



# Nucleus Manipulators



**Bechert**  
 Lens Nucleus Rotator  
 Blunt forked tip  
 Angled at 7mm  
 NM-2102



**Jaffe-Bechert**  
 Lens nucleus Rotator  
 Blunt forked tip  
 Angled at 9mm  
 NM-2103



**Clayman**  
 Nucleus Rotator  
 0.5mm wide with  
 V-notched tip  
 NM-2104



**Drysdale**  
 Nucleus Manipulator  
 Paddle shaped tip  
 NM-2105



**Pisacano**  
 Nucleus Rotator  
 Flat blunt tip  
 NM-2106



**Siepser**  
 Endocapsular Controller  
 Trumpet shaped tip  
 NM-2107



**Tennant**  
 Nucleus Manipulator  
 Ball shaped tip  
 NM-2108



# Spatulas & Loops

**Elsching**  
Cyclodialysis Spatula  
Flat, 1mm wide blade  
SL-2109



**Castroviejo**  
Synechiaie Spatula  
semi-sharp tips  
0.5mm wide  
double-end  
SL-2110




**Castroviejo**  
Cyclodialysis Spatula  
double-end  
0.25mm wide  
0.5mm wide  
SL-2111A  
SL-2111B



**Kimura**  
Platinum Spatula  
SL-2112



**Lewis**  
Lens Loop  
Small  
Large  
SL-2113S  
SL-2113L



**New Orleans**  
Lens Loop  
SL-2114



**Wilder**  
Lens Loop  
serrated  
SL-2115




# Punches & Polishers



## Gass

Corneoscleral punch  
1.5mm dia. Tip  
rotatable for  
Convenient positioning  
PP-2116



## Kelly

Descemet's  
Membrane Punch  
1.0mm dia.tip  
PP-2117



## Meyerhoefer

Chalazion Curette  
Size 0, 1.5mm dia.  
Size 1, 1.75mm dia.  
Size 2, 2.25mm dia.  
Size 3, 2.5mm dia.  
Size 4, 3.5mm dia.



PP-2118A PP-2118B PP-2118C PP-2118D

## Gills

Capsule Polishing Curette  
1.75mm dia. Cup  
with sandblasted base  
gentle backward angle  
PP-2119



## Nightingale

Capsule Polishing Curette  
1.75mm dia. Open ring  
with semi-sharp edge  
PP-2120



## Anis

Capsule Scrubber  
Gently curved flexible shank  
With roughened tip  
Ball tip  
Disc tip  
Reverse ball tip



PP-2121A PP-2121B PP-2121C

# Spatulas & Manipulators

**Lippman**  
Sideport Nucleus Spatula  
0.5mm dia. Round  
blunt forked tip  
**SM-2122**



**Girard**  
Synechiaie Spatula  
Semi-sharp tip  
0.5mm wide  
**SM-2123**



**Obstbaum**  
Lens Spatula  
Angled blade  
0.25mm wide  
**SM-2124**



**Hirschman**  
Lens Spatula  
Flat blade, 0.5mm wide  
With small notches for  
Manipulating IOLs  
**SM-2125**



**Jaffe**  
Lens Spatula  
Flat, smooth  
0.5mm wide  
**SM-2126**



**Lester**  
IOL Manipulator  
For Lenses with or  
Without hole  
Straight, with guard  
angled  
**SM-2127S**  
**SM-2127A**







**Kirby**  
 Lens Loop  
 And Muscle Hook  
 L -2128



**Gills-Welsh**  
 Lens Loop and Probe  
 with serrations  
 L -2129



**Morrison**  
 Lens Loop  
 And Probe  
 Very thin  
 L -2130



**Simcoe**  
 Lens Nucleus Loop  
 Fine serrations, strong curve  
 Double-end  
 L -2131



**Kansas**  
 Lens Loop, 3mm wide  
 Serrated, for small incision  
 Nucleus fragment  
 Removal  
 L -2132



**Keener-Arit**  
 Lens Loop  
 3mm wide, serrated  
 strong curve for  
 small incision nucleus  
 fragment removal  
 L -2133




**Keener**  
 Nucleus Divider  
 With two pre-formed  
 Snare wire loops  
 Complete  
 Wire loops(pk/3)  
 L -2134



# Spoons & Iris retractors

**Bunge**  
Evisceration Spoon  
Small  
Large  
SI-2134S  
SI-2134L




**Wells**  
Enucleation Spoon  
SI-2135




**Rizzuti**  
Iris Retractor  
And Lens  
Expressor  
2mm wide lip  
SI-2136



**Azar**  
Iris Retractor  
Autoclavable plastic  
Left  
Right  
SI-2137L  
SI-2137R



**Graether**  
Collar Button  
Micro iris retractor  
Straight, with guard  
angled  
SI-2138S  
SI-2138A



# Nucleus Dividers & Spoons

**Kansas**  
**Nucleus Trisector**  
 For manual small incision surgery  
 For right handed  
 Surgeons  
 For left handed  
 Surgeons  
**NS-2139**



**Kansas**  
**Nucleus Bisector**  
 For manual small incision surgery  
**NS-2140**



**Kansas**  
**Nucleus Vectis, solid**  
 Used as a "cutting block"  
 For Kansas nucleus  
 Bisector and trisector  
 For right handed surgeons  
 For left handed surgeons  
**NS-2141**



**Culler**  
 Lens Spoon  
 Sterling Silver  
**NS-2142**



**Barraquer**  
 Lenticular Spoon  
 For corneal button transfer  
**NS-2143**



**Paton**  
 Spatula and Spoon  
 For corneal button transfer  
**NS-2144**



# Iris Hooks

**Graether**  
 Collar Button  
 Micro iris retractor  
 And posterior Capsule polisher  
 With side irrigation  
 Straight  
 Angled  
 IH-2145S  
 IH-2145A



**Bonn**  
 Micro Iris Hook  
 Blunt, round handle  
 With protective Guard  
 IH-2146




**Bonn**  
 Micro Iris Hook  
 Blunt, flat handle  
 IH-2147



**Jaffe-Knolle**  
 Iris Hook, very blunt  
 Small, angled  
 IH-2148



**Knolle**  
 Iris Hook, very blunt  
 Small, straight  
 IH-2149



**Shepard**  
 Iris Hook  
 Small, blunt  
 IH-2150






### Sinskey

Lens Manipulating Hook  
With round handle  
0.25mm dia. Blunt tip  
straight, with guard  
angled

IH-2151S  
IH-2151A



### Sinskey II

Lens manipulating Hook  
With round handle  
0.2mm dia. Blunt tip  
straight  
angled


IH-2152S  
IH-2152A



### Sinskey II

Lens Manipulating Hook  
With round handle  
0.2mm dia. Blunt tip  
straight, with guard  
angled

IH-2153S  
IH-2153A



### Lewcky

Lens manipulating Hook  
0.15mm dia. Blunt tip  
angled

IH-2154



### Hirschman

Iris Hook  
Blunt, angled  
IH-2155



# Hooks

## Maltzman-Fenzl

Lens Manipulating Hook  
0.15mm dia. V-shaped tip  
round handle  
straight, with guard  
angled  
H -2156



## Galand

Lens Manipulating Hook  
angled  
H -2157



## Jaffe-Maltzman

Lens Manipulating Hook  
0.15mm dia. V-shaped tip  
angled  
H -2158



## Osher

Y-Hook  
for manipulating IOLs  
straight, with guard  
H -2159



## Fixation Forceps

double, sharp  
small  
large  
H-2160S  
H-2160L



## Twist Hook

for scleral fixation  
right  
left  
H-2161R  
H-2161L



## Anis

Scleral Twist Fixation Cone  
2.5mm dia. With three  
teeth, with guard  
H -2162



# Hooks



**Kuglen**  
Iris Hook  
And Lens Manipulator  
Push-pull model  
Straight  
Angled  
H-2163S  
H-2163A



**Kuglen**  
Iris Hook  
And Lens Manipulator  
Clover leaf model  
Straight  
H -2164



**Kuglen**  
Iris Hook  
And Lens manipulator  
Push-pull model  
With round handle  
Straight, with guard  
Angled  
H-2165S  
H-2165A



**Rentsch"Boat Hook"**  
Lens Manipulator  
This multi-purpose  
Instrument may be used as  
An iris/capsule retractor or  
Pusher/manipulator for IOLs  
With and without holes  
Straight  
With guard  
H -2166



**Fenzi**  
Lens Manipulating Hook  
0.15mm dia. T-shaped tip  
straight,  
with guard  
angled  
H-2167S  
H-2167A



# Needle Holders

## McPherson

Needle Holder

very delicate jaws

straight, with lock

straight, without lock

curved, with lock

curved, without lock

NH-2168A NH-2168C

NH-2168B NH-2168D



## Jaffe

Needle Holder

delicate conical jaws

straight, with lock

curved, with lock

curved, without lock

NH-2169A

NH-2169B

NH-2169C



## Castroviejo

Needle Holder

delicate jaws

straight, with lock

straight, without lock

curved, with lock

curved, without lock

NH-2170A

NH-2170B

NH-2170C

NH-2170D







# Needle Holders



**Castroviejo**  
Needle Holder  
standard jaws, with lock  
straight  
curved  
NH-2171A  
NH-2171B  
NH-2171C



**Kalt**  
Needle Holder  
straight smooth jaws  
delicate  
NH-2172



**Drews**  
Needle Holder  
extra delicate jaws  
straight, with lock  
straight, without lock  
curved, with lock  
curved, without lock  
NH-2173A  
NH-2173B  
NH-2173C  
NH-2173D



**Barraquer**  
Needle Holder  
short model  
very delicate jaws  
straight, with lock  
straight, without lock  
curved, with lock  
curved, without lock  
NH-2174A NH-2174C  
NH-2174B NH-2174D



# Needle Holders


## Anis

Needle Holder  
 extra delicate jaws  
 straight, with lock  
 straight, without lock  
 curved, with lock  
 curved, without lock  
**NH-2175S**  
**NH-2175C**



## Katena

Needle Holder  
 extra delicate jaws  
 straight, with lock  
 straight, without lock  
 curved, with lock  
 curved, without lock  
**NH-2176A**  
**NH-2176B**  
**NH-2176C**  
**NH-2176D**




## McInyre

Fish Hook Needle Holder  
 extra delicate jaws  
 straight, with lock  
**NH-2177**



# Needle Holders

**Troutman**  
Needle Holder  
delicate, strongly curved jaws,  
with lock  
without lock  
**NH-2178A**  
**NH-2178B**




**Barraquer**  
Needle Holder  
very delicate curved jaws  
with lock  
without lock  
**NH-2179A**  
**NH-2179B**



**Barraquer**  
Needle Holder  
delicate jaws  
straight, with lock  
straight, without lock  
curved, with lock  
curved, without lock  
**NH-2180A**  
**NH-2180B**  
**NH-2180C**  
**NH-2180D**



**Barraquer**  
Needle Holder  
short model  
standard jaws, curved  
with lock  
without lock  
**NH-2181**



# Speculums



**Barraquer**  
Wire Speculum  
small  
large, heavy wire  
SP-2182S  
SP-2182L



**Barraquer**  
Wire Speculum  
with solid blades  
small  
large  
SP-2183S  
SP-2183L




**Grandon-Barraquer**  
Wire Speculum  
with suture notches or holes  
to obtain additional traction  
for the upper lid  
wire blades  
solid blades  
SP-2184W  
SP-2184S



**Feaster-Barraquer**  
Wire Speculum  
wide, solid blades  
SP-2185



**McIntyer**  
Wire Speculum with  
gently curved spring.  
V-shaped blades help to retain  
disposable drape under lids  
close wire  
open wire  
SP-2186




# Speculums



**Kartz-Barraquer**  
Wire Speculum  
with open blades  
for unobstruted access  
to the cataract incision  
with phaco or I/A tip  
heavy wire, large  
standard wire, large  
standard wire, small  
**SP-2187L**  
**SP-2187S**



**Kartz-Barraquer**  
Wire Speculum with wings  
for retracting surgical  
drape while allowing  
unobstruted access  
to the cataract incision  
with phaco or I/A tip  
standard wire, small  
**SP-2188W**  
**SP-2188S**



**Tennant**  
Lid Retractor  
may be clamped to surgical  
drape with suture or elastic band  
pair  
**SP-2189**



**Sauer**  
Lid Speculum  
for infants  
for premature infants  
**SP-2190**



**Infant**  
wire Speculum  
with gently curved springs  
4mm wide blades  
**SP-2191**

**Cook**  
Eye Speculum with locking screw  
solid,small #1  
solid,medium #2  
solid,large #3  
fenestrated, medium #2  
fenestrated, large #3  
**SP-2192S**  
**SP-2192M SP-2192L**




# Speculums


**Castroviejo**  
 Eye Speculum  
 with finely controlled opening  
 and closing mechanism  
 small #1  
 SP-2193




**Castroviejo**  
 Eye Speculum  
 with finely controlled opening  
 and closing mechanism  
 medium #2  
 SP-2194



**Fishkind-castroviejo**  
 Lid Speculum with open blades  
 for use with phaco I/A tip  
 SP-2195



**McPherson**  
 Eye Speculum  
 solid blades with  
 vertical suture posts  
 SP-2196




# Speculums

**Williams  
Eye Speculum  
open blade, small  
SP-2197**



**Williams  
Eye Speculum  
open blade, large  
SP-2198**



**Lancaster  
Eye Speculum  
solid blade  
SP-2199**




**Lester-Burcs  
Eye Speculum  
solid blade with  
flanges for retaining  
eye lashes  
SP-2200**



# Speculums


## Meaumennee-park

Eye Speculum  
with canthus hook,  
blade locking screws may  
be positioned vertically  
or horizontally  
solid blades  
fenestrated blades  
SP-2201




## Kuglen

Eye Speculum  
with canthus hook,  
solid blades with  
swivel ball joints  
and locking screws  
SP-2202



## Guyton-park

Eye Speculum  
with vertical suture posts  
and canthus hook  
solid blades  
fenestrated blades  
SP-2203





# Speculums


## Smirmaul

Eye Speculum  
 individually adjustable  
 blade for elevating lids  
 solid blades  
 wire blades  
 SP-2204




## Schott

Eye Speculum  
 with rotatable adjustable  
 for precise control of  
 blades elevation  
 wire blades  
 SP-2205



## Flieringe

Scleral Fixation Ring  
 set of 8 sizes (15-22mm)  
 14mm dia.  
 15mm dia.  
 16mm dia.  
 17mm dia.  
 18mm dia.  
 19mm dia.  
 20mm dia.  
 21mm dia.  
 22mm dia.  
 SP-2206



# Retractors

**Jaeger**  
Lid plate  
stainless steel  
RE-2207



**Wright**  
Fascia Needle  
with 1mm x 6mm  
oval hole  
14cm long  
RE-2208



**Knapp**  
Lacrimal Retractor  
4-prong  
blunt  
sharp  
RE-2209



**Stevenson**  
Lacrimal Sac Retractor  
3 x 3 blunt prongs  
self-retaining  
adjustable spread  
RE-2210



**Stevenson**  
Lacrimal Sac Retractor  
solid blades with  
serrated edge  
self-retaining  
adjustable spread  
RE-2211



**Agricola**  
Lacrimal Sac Retractor  
3 x 3 sharp prongs  
RE-2212



# Retractors



**Schepens**  
Froked Orbital Retractor  
SP-2213




**Arruga**  
Globe Retractor  
SP-2214



**Jaffe**  
Wire Lid Retractor  
may be clamped to  
surgical drape with  
suture or elastic band  
pair, large  
pair, small  
SP-2215



**Desmarres**  
Lid Retractor  
solid blade  
size 0  
size 1  
size 2  
size 3  
SP-2216A  
SP-2216B  
SP-2216C  
SP-2216D



**Jaeger**  
Lid Plate  
clear plastic  
SP-2217



# Knives



**Ziegler**  
 Knife Needle  
 3mm  
 5mm  
 KN-2218



**Girard-swain**  
 Knife Needle  
 1.5mm cutting edge  
 5mm  
 KN-2219



**Wheeler**  
 Discission Knife  
 18mm cutting edge  
 KN-2220



**Emery**  
 cystotome  
 double-edge  
 for push and pull cutting  
 blade, straight  
 blade, angled  
 chuck handle  
 KN-2221



**Stiletto Blade**  
 For micro-vitro-retinal  
 stab incisions, lancet  
 shaped, 19 gauge  
 1.2mm wide  
 blade  
 chuck handle  
 KN-2222



**Paufique**  
 Graft Knife  
 angled tip  
 curved cutting edge  
 KN-2223



# Knives & Spuds

## Tooke

Coeneal Knife  
straight blade with  
curved cutting edge  
KS-2224



## Gill

Coeneal Knife  
curved cutting edge  
straight  
curved  
KS-2225



## Golf-Club

Foreign Body Spud  
angled, lightly rounded tip  
KS-2226



## Francis

Foreign Body Spud  
angled, flag shaped tip  
KS-2227



## Ellis

Foreign Body Spud  
very thin, lightly  
curved and rounded tip  
KS-2228



## Spud & Gouge

In reversible handle  
for removing foreign bodies  
KS-2229




## Wills Eye Hospital

Foreign body spud & bar  
delicate curved spud on  
one end and small barr  
for rust ring reoval  
on opposite end  
KS-2230



# Cannulas



**West**  
Irrigating Cannula  
Blunt tip, 23 gauge  
Side Opening  
KS-2231



**Kelman**  
Cyclodialysis cannula  
spatulated tip  
side opening, 19 gauge  
left opening  
Right opening  
KS-2232



**Air injection**  
Cannula, angled 5mm from tip  
30 gauge  
27 gauge  
KS-2233




**Sheets**  
Cyclodialysis Cannula  
spatulated tip  
one front and two side  
openings, 19 gauge  
KS-2234



**Viscoelastic**  
Aspirating cannula  
Angled 10mm from tip  
22 gauge  
19 gauge  
KS-2235



**Troutman**  
Alpha-Chymotrypsin Cannula  
delicate olive tip  
25 gauge  
KS-2236



**Rainin**  
Air injection cannula  
Spatulated tip  
Angled at 7mm  
30 gauge  
27 gauge  
KS-2237



**Wfesh**  
Flat Olive tip Cannula  
25 gauge  
straight  
left  
right  
KS-2239S  
KS-2239L  
KS-2239R




**Randolph**  
Cyclodialysis cannula  
Curved spatulated tip  
Angled at 13mm  
21 gauge  
19 gauge  
KS-2238



# Blade Breakers



**Aldgerbrush**  
 Corneal Rust  
 Ring Remover  
 battery powered  
 complete  
 extra burr 0.5mm  
 extra burr 1.0mm  
**KS-2240**



**Castroviejo**  
 Blade Breaker  
 and holder  
 with sliding lock  
**KS-2241**



**Castroviejo**  
 Blade Breaker and holder  
 short model  
 pointed jaws  
 squared jaws  
**KS-2242**



**Swiss Model**  
 Blade Breaker  
 and holder  
 spherical tip  
**KS-2243**




**Katena**  
 Double Edge Blades  
 for blade breakers  
 pk/10 non-sterile  
 pk/50 sterile  
**KS-2244**



**Elliot**  
 Scleral Trephine  
 1.0mm dia.  
 1.5mm dia.  
 2.0mm dia.  
 Handel  
**KS-2245**




# I-A Cannulas



## Simcoe

Irrigating-Aspirating Cannula  
Original model, with aspiration  
Through tubing hub and irrigation  
Through luer-lock hub. 23 gauge  
Thin wall, 15mm long curved  
Tubes, with 0.3mm dia. Aspirating  
Port facing into curve and slightly  
Shorter irrigating tip with 45°  
Angled front opening. Designed  
For true Simcoe I/A technique.


IA-2246



## Swiss

Irrigating-Aspirating cannula  
for aspiration through tubing hub  
and irrigation through luer-lock  
hub. 23 gauge thin wall, 15mm long  
curved tubes, 0.3mm dia. Port on left  
side of aspirating tube. Slightly  
angled front opening.


IA-2247



## Reversed simcoe

Irrigating-Aspiratin Cannula  
Same as K7-4300 but tubes reversed.  
For surgeons who like the simcoe I/A  
Cannula, but prefer to use a more  
conventional technique of aspirating  
through the luer-lock hub cannula .

IA-2248



## Ambrosio

Irrigating-Aspirating  
Cannula For aspirating  
through tubing hub and  
Irrigation through luer-lock  
hub. 23 gauge thin wall,  
15mm long, curved,  
front openings.

smooth tip  
rough tip for polishing

IA-2249





# Cystotomes




**Kelman**  
Cystotome  
With sharp point and  
Sharp cutting edge  
22 gauge  
IA-2250



**Knolle-Kelman**  
Cystotome  
small blade  
with sharp point and  
sharp cutting edge  
22 gauge  
IA-2251



**Kratz**  
Cystotome  
With sharp point and  
Sharp cutting edge  
Vaulted, 22 gauge  
IA-2252



**McIntyre**  
Cystotome  
angled 10mm from tip  
with sharp point, 23 gauge  
IA-2253



**Kelman**  
Cystotome  
With sharp point and  
Blunt blade, 22 gauge  
IA-2254



**Lieppman**  
Micro Cystotome  
extra small, with sharp point  
and sharp cutting edge  
22 gauge  
IA-2255



**Kelman**  
Double Blade Cystotome  
With sharp points and  
Blunt blades, 22 gauge  
IA-2256



**Graether**  
Collar Button  
micro iris retractir  
and posterior capsule polisher  
with side irrigation  
straight  
angled  
IA-2257



**Double-Edge**  
Cystotome  
With sharp point and  
sharp cutting edges  
22 gauge  
IA-2258



**Shammas**  
Posterior Capsule  
Polishing Curette  
1mm diameter cup with  
sandblasted base  
22 gauge  
IA-2259

# Scalpel Handle & Blades



**BH-2274A Miniature Blade**  
**BH-2274B** Sterile, disposable  
**BH-2274C** #63 (pk/6)  
**BH-2274D** #64 (pk/12)  
**BH-2274E** #65 (pk/12)  
**BH-2274F** #67 (pk/12)  
 #69 (pk/6)  
 #69B (pk/6)



**Lamellar Blade**  
 sterile, disposable (pk/3)  
 #66-60  
 #66F  
**BH-2273A**  
**BH-2273B**



**Super Blade**  
 sterile, disposable (pk/12)  
 #76-15,15\*  
 #76-30,30\*  
 #76-45,45\*  
**BH-2275A**  
**BH-2275B**  
**BH-2275C**

**Chuck Handle**  
 #3H hexagon 10cm  
 #3K round 10cm  
 #3KA round 7.5cm  
 #3KB round 5.0cm

**BH-2276A**  
**BH-2276B**  
**BH-2276C**  
**BH-2276D**



## Stiletto Knife

For Micro-Vitreo-Retinal  
 Stab incisions  
 Sterile, disposable (pk/3)  
 20 gauge  
 19 gauge  
**BH-2277**



## Bard-Parker Handle

For disposable slot blades  
 #3 **BH-2278A**  
 #9 **BH-2278B**



## Bard-Parker Blade

Sterile, disposable (pk/50)  
 #11 **BH-2279A**  
 #15 **BH-2279B**




# Capsule Polishers



**Drews**  
Posterior Capsule Polisher  
Lightly angled, 2mm I.D loop  
With sandblasted base  
0.2mm irrigating port facing  
into circle, 25 gauge  
IA-2263



**Jensen**  
Posterior Capsule Polisher  
carbide impregnated on  
Posterior surface of olive tip  
front opening, gently curved  
25 gauge  
23 gauge  
IA-2264




**Kraff**  
Posterior Capsule  
Polishing Curette  
1.75mm diameter cup  
with semi-sharp edge.  
Irrigating port in center  
Of cup, 25 gauge  
IA-2265



**Jensen**  
Capsule Polisher  
carbide impregnated  
All around olive tip  
front opening, gently curved  
25 gauge  
23 gauge  
IA-2266



**Kratz**  
Scratcher  
Front opening with  
Carbide impregnated tip  
Lightly angled, 22 gauge,  
IA-2267



**Knolle**  
Posterior Capsule Polisher  
2mm angled tip, front opening  
finely sandblasted  
30 gauge  
23 gauge  
IA-2268



**Rubman-katzin**  
Posterior Capsule Polisher  
Carbide impregnated tip  
With smooth closed end.  
Irrigating side port facing  
Into curve, 22 gauge  
IA-2269



**Simcoe**  
Cortex Extractor and Capsule Polisher  
side opening facing into curve  
23 gauge, sandblasted  
13mm  
21mm  
IA-2270




**Shapiro**  
Posterior Capsule Polisher  
Carbide impregnated tip  
With dual 0.3mm sideports  
curved, 23 gauge  
IA-2271




**Hydro-Dissection**  
Cannula, 25 gauge  
Angled flat tip  
IA-2272

# Disposable Cannulas



**Bishop-Harmon**  
Irrigating Cannula  
Flat tip, angled at  
8mm from the tip  
sterile, disposable  
19 gauge  
DC-2280



**Hydrodissection**  
Cannula, flat tip  
Angled at 7mm  
From the tip  
Sterile, disposable  
25 gauge  
DC-2281




**Hydrodissection**  
cannula, U-Shaped  
with flat tip  
sterile, disposable  
25 gauge  
DC-2283



**J-Shape Cannula**  
Sterile, disposable  
Straight, 25 gauge  
angled, 25 gauge  
straight, 23 gauge  
angled, 26 gauge  
DC-2285




**Simcoe**  
Cortex Extractor  
23 gauge X 13mm  
0.3mm aspiration port  
sterile, disposable  
DC-2287



**Simcoe**  
Irrigating-Aspiration  
Cannula, 23 gauge thin wall  
0.3mm aspiration port  
with tubing and adaptor  
sterile, disposable  
DC-2282



**Kartz**  
Capsule Polisher  
end opening, angled at  
7mm from the tip  
sterile, disposable  
23 gauge  
21 gauge  
DC-2284



**Kartz**  
Capsule Polisher  
close tip with side port, angled  
at 7mm from the tip  
sterile, disposable  
23 gauge  
21 gauge  
DC-2286



**Silicone**  
Tip Cannula  
angled at 7mm from  
the tip, 27 gauge  
sterile, disposable  
DC-2288



# Bipolar Forceps




**Bipolar Forceps**  
 Jeweler's Type  
 pointed tips, straight  
 plain  
 insulated  
 19 gauge  
 DI-2289



**Bipolar Forceps**  
 McPherson Type  
 lightly blunted tips  
 straight  
 DI-2290




**Bipolar Forceps**  
 McPherson Type  
 lightly blunted tips  
 angled  
 DI-2291




**Bipolar**  
 Coaptation Forceps  
 pointed tips  
 DI-2292



**Bipolar Cord**  
 disposable(10/box)  
 reusable  
 DI-2293




# Cannulas




**McIntyre-Binkhorst**  
Irrigating Cannula  
smooth blunt tip  
front opening, 26 gauge  
right  
straight  
left  
C-2260R  
C-2260S  
C-2260L




**Gimbel**  
Hydro-dissection Cannula  
smooth blunt tip  
front opening, 30 gauge  
C -2263



**Knolle-Pearce**  
Irrigating Vectis  
dual fluid flow through  
three front opening ports  
23 gauge loop  
C -2264



**Drews**  
Irrigating Vectis  
dual fluid flow through  
one front and two backward  
opening ports, 23 gauge loop  
C -2266



**Sheets**  
Irrigating vectis  
dual fluid flow through  
three large front opening  
ports, serrated 20 gauge loop  
C -2261



**Hofmann**  
Nucleus loops  
3mm wide, serrated  
extra delicate  
non-irrigating  
C -2262



**Simcoe**  
Irrigating Lens Nucleus Loop  
two front opening ports  
finely serrated, strong curve  
25 gauge  
right  
left  
C-2265R  
C-2265L



**Anis**  
Irrigating Vectis  
dual fluid flow through  
two front opening ports  
23 gauges heart shaped loop  
enables surgeon to retract  
iris and capsular flap  
while irrigating  
C -2267



# Cannulas



**Nishi**  
Nucleus Delivery Cannula  
flat insertion plate with  
central irrigating port.  
Restraining fin limits  
depth of insertion and  
may be used to apply gentle  
pressure to the sclera.  
Tip of insertion plate is  
placed under the nucleus,  
while continuous irrigation  
helps to float the nucleus  
forward across the plate,  
reducing zonular stress.  
Right  
left  
C -2268



**Saiz**  
Irrigating Cassule Retractor  
20 gauge, straight  
C - 2269



**Kratz**  
Lens Manipulating Cannula  
0.2mm diameter tip  
frot opening, 30 gauge  
C -2270



**Thurmond**  
Nucleus Irrigating Cannula  
and Iris Retractor, front opening  
smooth downward curved tip  
19 gauge  
C -2271



**Osher**  
Lens Manipulating Cannula  
with sideport and pigtail tip  
24 gauge  
C -2272




**Nishi**  
Irrigating Capsule Rectors  
for hydro-dissection of the  
lens nucleus  
C -2273



**Kuglen**  
Irrigating Lens Manipulator  
disc shaped tip with 0.25mm  
diameter peg protruding  
from its center, 25 gauge  
straight  
angled  
C-2274S  
C-2274A

# Cataract Set



## PLANNED EXTRA-CAPSULARA CATARACT EXTRACTION

### CS-2275

- |                                 |  |
|---------------------------------|--|
| BSS Irrigator                   | Utrata Capsulorhexis Forceps           |
| Towel Clamp                     | IoL Forceps                            |
| Mosquito Forceps                | IoL Manipulator                        |
| Serrefine (2)                   | Kuglen Iris Hook                       |
| Wire Speculum                   | Cyclodialysis Spatula                  |
| Super Sharp Blade               | Needle Holder, micro                   |
| Blade Handle                    | Needle Holder, stright                 |
| Corneal Section Scissors, right | Lens Loops & Probe                     |
| Corneal Section Scissors, Left  | Thornton Limbal Ruler                  |
| Westcott Scissors, blunt        | Kansas Nucleus Trisector               |
| Stitch Scissors                 | Kansas Nucleus Bisecto                 |
| Capsulotomy Scissors            | Kansas Vectis, solid                   |
| Iris Scissors                   | Kansas Nucleus Loop                    |
| Caliper                         | Kansas Nucleus Divider                 |
| Conjunctive Forceps             | Loop for Nucleus Divider               |
| Superior Rectus Forceps         | Kansas Nucleus Fragment Forceps        |
| Colibri Forceps 0.12mm          | Viscoelastic Aspirating Cannula, 22ga. |
| Diamond Step Knife              | Air Injection Cannula, 30ga.           |
| Intraoperative Keratometer      | Viscoelastic Injection Cannula, 26ga.  |
| Tissue Forceps                  | Cystotome                              |
| Suturing Forceps, 0.12mm        | Capsule Polisher                       |
| Tying Forceps, Stright          | I/A Cannula                            |
| Tying Forceps, angled           | Sterilizing Case(not illustrated)      |
| Kelman-mcPherson Forceps        | K-Spong Spears(not illustrated)        |






# Enucleation Set



**ENUCLEATION**  
**Speculum**  
**Tissue Forceps**  
**Mosquito Forceps**  
**Tenotomy Scissors**  
**Enucleation Scissors**  
**Muscle Hook**  
**Sterilizing Case(not illustrated)**  
**ES-2276**



# Radial Keratotomy Set




## RADIAL KERATOTOMY

### RKS-2276

- Wire Speculum, solid
- Diamond Micrometer, double edge
- Blade Gauge
- Marking Pen
- Thornton Fixation Ring
- Fixation Forceps, straight
- Fixation Forceps, angled
- Visual Axis Marker
- Optic Zone Maker, 3.0mm
- Optic Zone Maker, 3.25mm
- Optic Zone Maker, 3.5mm
- Optic Zone Maker, 3.75mm

- Optic Zone Maker, 4.0mm
- Optic Zone Maker, 4.25mm
- Optic Zone Maker, 4.5mm
- Optic Zone Maker, 5.0mm
- Radial Marker 6 Line
- Radial Marker 8 Line
- Incision Spreading Forceps
- Incision Depth Gauge
- Irrigating Cannula, 27ga
- Sterilizing Case (not illustrated)
- K-Spong Spears (not illustrated)





## LID SURGERY


### LS-2277

BSS Irrigator  
 Towel Clamp(2)  
 Mosquito Forceps, curved(2)  
 Tenotomy Scissors  
 Utility Scissors  
 Miniature Blade#69  
 Blade Handle  
 Westcott Scissors, blunt  
 Stitch Scissors  
 Ptosis Forceps  
 Needle Holder, delicate  
 Kalt Needle Holde  
 Lid Plate  
 Caliper  
 Lid Retractor #1  
 Lid Retractor #2  
 Lid Retractor #3

Muscle Hook#1  
 Muscle Hook#1  
 Fixation Hook, sharp  
 Chalazion Curette #1  
 Chalazion Curette #2  
 Chalazion Curette #3  
 Chalazion Curette #4  
 Fixation Frceps  
 Conjunctive Forceps  
 Tissue Forceps  
 Fixation Forceps,angled  
 Suturing Forceps  
 Tying Forceps, curved  
 Chalazion Forceps, round  
 Chalazion Forceps,ovel  
 Lid Forceps  
 Sterilizing Case(not illustrated)  
 K-Spong Spears(not illustrated)



# Lacrimal Set




## LACRIMAL SURGERY

### LS-2278 (Dacryocystorhinostomy & Intubation)

- |                               |                                 |
|-------------------------------|---------------------------------|
| BSS Irrigator                 | Lacrimal Chisel, straight       |
| Towel Clamp(2)                | Pigtail Probe                   |
| Mosquito Forceps, straight(2) | Retrieving Hook                 |
| Mosquito Forceps, curved(2)   | Intubation set                  |
| Kerrison Rongeur              | Lacrimal probe set              |
| Nasal Speculum                | Dressing Forceps                |
| Bone Rongeur                  | Tissue Forceps                  |
| Mallet                        | Utility Forceps                 |
| Maniature Blade #65           | suturing Forceps                |
| Blade Handle                  | Rake Retractor(2)               |
| Self Retaining Retractor      | Stevns Scissors                 |
| Lacrimal Cannula, straight    | Utility Scissors                |
| Lacrimal Cannula, curved      | Needle Holder                   |
| Lacrimal Dilator              | Stitch Scissors                 |
| Periosteal Elevator           | K-Spong Spears(not illustrated) |
| Lacrimal Chisel, straight     |                                 |



# Enucleation Set



## ASTIGMATIC KERATOTOMY

### ES-2279

- Wire Speculum, solid
- Diamond Micrometer, double edge
- Blade Gauge
- Marking Pen
- Thornton Fixation Ring
- Fixation Forceps, double
- Fixation Forceps, 0.12mm
- Visual Axis Marker
- Degree Gauge

- Aix Marker
- T-Incision Marker
- Optic Zone Maker
- Thornton Corneal Ruler
- Arcuate Incision Marker
- Irrigating Cannula, 27 ga.
- Incision Spreading Forceps
- Incision Depth Gauge
- Irrigating Cannula, 27ga
- Sterilizing Case (not illustrated)
- K-Spong Spears (not illustrated)



# Sterilizing Case




**Silicone Sterilizing Case**  
Venting Holes  
26cm x 15cm x 2.5cm  
**SC-2280**



**Silicone Sterilizing Case**  
Venting Holes  
26cm x 24cm x 2.5cm  
**SC-2281**



**Silicone Sterilizing Case**  
Venting Holes  
Full silicone finger mat top & bottom  
26cm x 24cm x 2.5cm  
accomadetes 16 instruments  
**SC-2282**




**Silicone Sterilizing Case**  
Venting Holes  
26cm x 15cm x 2.5cm  
accomadetes 10 instruments  
**SC-2283**







# Titanium Tissue Forceps



ADSON  
BD-30T  
120 mm




ADSON  
BD-31T  
120 mm




ADSON-BROWN  
BD-141T  
120 mm


9x9 Teeth




BD-09T  
130 mm




BD-75T  
145 mm



BD-16T  
160 mm




COLLEGE  
TT-01  
150 mm




MERIAM  
TT-02  
160 mm


# Titanium Tissue Forceps




GERALD  
MF-101T  
180 mm



GERALD  
MF-102T  
180 mm



GERALD  
MF-103T  
180 mm




GERALD  
MF-104T  
180 mm






# Titanium Tissue Forceps




BD-32T BD-33T

BD-32T/BD-33T  
100 mm




BD-35T BD-36T

BD-35T/BD-36T  
100 mm




BD-60T BD-62T

BD-60T/BD-62T  
110 mm




BD-42T BD-43T


BD-42T/BD-43T  
110 mm




BD-285T  
90 mm



BD-270T  
130 mm



BD-341T  
115 mm




SEMKEN  
BD-46T  
130 mm

1x2 Teeth

# Titanium Tissue Forceps

Mikro-Pinzette  
Ø 8 mm



0,4 mm  
1/1  
**MF-130T**  
150 mm

0,4 mm  
1/1  
**MF-131T**  
150 mm

0,8 mm  
1/1  
**MF-132T**  
150 mm

0,8 mm  
1/1  
**MF-133T**  
150 mm

0,4 mm  
1/1  
**MF-134T**  
180 mm

0,4 mm  
1/1  
**MF-135T**  
180 mm

0,8 mm  
1/1  
**MF-136T**  
180 mm

0,8 mm  
1/1  
**MF-137T**  
180 mm

0,4 mm  
1/1  
**MF-138T**  
210 mm


0,4 mm  
1/1  
**MF-139T**  
210 mm

0,8 mm  
1/1  
**MF-140T**  
210 mm

0,8 mm  
1/1  
**MF-141T**  
210 mm

MF-130T/MF-141T

Mikro-Knüpfpinzette  
Ø 8 mm



6 x 0,4 mm  
1/1  
**MF-362T**  
150 mm

6 x 0,4 mm  
1/1  
**MF-363T**  
150 mm

6 x 0,8 mm  
1/1  
**MF-364T**  
180 mm

6 x 0,8 mm  
1/1  
**MF-365T**  
180 mm

6 x 0,8 mm  
1/1  
**MF-366T**  
210 mm

6 x 0,8 mm  
1/1  
**MF-367T**  
210 mm


MF-362T/MF-367T



# Titanium Ring Forceps



Mikro-Ring Forceps  
Ø 8 mm



MF-142T/MF-347T

Smooth

Serrated

1 x 0,5 mm

2 x 1 mm

1 x 0,5 mm

2 x 1 mm



MF-142T  
150 mm

MF-145T  
150 mm

MF-342T  
150 mm

MF-345T  
150 mm

1 x 0,5 mm

2 x 1 mm

1 x 0,5 mm

2 x 1 mm



MF-143T  
180 mm

MF-146T  
180 mm

MF-343T  
180 mm

MF-346T  
180 mm

1 x 0,5 mm

2 x 1 mm

1 x 0,5 mm

2 x 1 mm



MF-144T  
210 mm


MF-147T  
210 mm

MF-344T  
210 mm

MF-347T  
210 mm

# Titanium Tissue Forceps

Atraumatic Micro Forcep



MF-348T/MF-359T

Atrauma-Fino  
1,2 mm



1/1

**MF-348T**  
150 mm

Atrauma-Fino  
1,2 mm



1/1

**MF-349T**  
150 mm

Atrauma-Standard  
1,8 mm



1/1

**MF-354T**  
150 mm

Atrauma-Standard  
1,8 mm



1/1

**MF-355T**  
150 mm

Atrauma-Fino  
1,2 mm



1/1

**MF-350T**  
180 mm

Atrauma-Fino  
1,2 mm



1/1

**MF-351T**  
180 mm

Atrauma-Standard  
1,8 mm



1/1

**MF-356T**  
180 mm

Atrauma-Standard  
1,8 mm



1/1

**MF-357T**  
180 mm

Atrauma-Fino  
1,2 mm



1/1

**MF-352T**  
210 mm

Atrauma-Fino  
1,2 mm



1/1

**MF-353T**  
210 mm

Atrauma-Standard  
1,8 mm



1/1

**MF-358T**  
210 mm

Atrauma-Standard  
1,8 mm




1/1


**MF-359T**  
210 mm




# Titanium Tissue Forceps




DE BAKEY  
BD-80T/BD-163T




BD-160T  
160 mm




BD-161T  
180 mm




BD-162T  
200 mm




BD-163T  
240 mm




BD-80T  
160 mm




BD-81T  
180 mm



BD-84T  
200 mm



BD-86T  
240 mm




DE BAKEY  
BD-82T  
160 mm


# Titanium Micro Scissors & Needle Holders

Micro Needle Holders


Micro Scissors



1/1 1/1  
FD-562T FD-563T



1/1 1/1  
FD-574T FD-575T



1/1 1/1  
FD-587T FD-588T


FD-562T/FD-563T  
180 mm

FD-574T/FD-575T  
180 mm


FD-587T/FD-588T  
180 mm




# Fine Artery Forceps




HALSTED-MOSQUITO  
P-401/P-402  
120 mm



HALSTED-MOSQUITO  
P-403/P-404  
120 mm



KELLY  
P-405/P-406  
140 mm



CRILE  
P-407/P-408  
140 mm

# Ceramic Coated Instruments



Castroveijo Caliper  
0042  
20mm



0043  
Castroveijo Scissor  
10cm CVD



0047  
Trout Barraquer  
Needle Holder  
13cm CVD



0048  
Vannas Scissor  
8cm CVD



0062-  
Castroveijo Scissor  
CVD 10cm



0069  
Westcott Tenotomy Scissor  
10cm CVD







# Ceramic Coated Forceps



0049  
McPherson Time Forcep  
10cm



0053  
MicroBishop Tissue Forcep  
8.5cm



0056  
Utrata Capsular Hexis Forcep  
Moria 11cm




0057  
Suture time Forcep  
Straight 7mm



# Ceramic Coated Instruments



0060  
Suture Time forcep  
CVD 7mm



0064  
Castroviejo Suture Forcep  
Straight 0.12mm



0065  
Castroviejo Suture Forcep  
Straight 0.12mm



0067  
Utrata Capsular Hexis Forcep  
11cm



0070- IOL Hook



0066  
Barraquer Eye Speculum  
4cm



## A

Arrowsmith.....	03
Anis.....	13, 16, 21, 28, 32, 52
Azar.....	24
Agricola.....	40
Arruga.....	41
Air Injection.....	44
Aldgerbrush.....	45
Ambrosio.....	46
Air Injection Cannula, 30ga.....	54
Arcuate Incision Marker.....	59

## B

Botvin.....	01
Bross.....	02
Bishop-Harmom.....	05
Barken.....	04
Barraquer.....	07, 18, 25, 31, 34, 32, 33
Biard.....	08
Backhaus.....	10
Baby jone.....	10
Berke.....	11
Bechert.....	19
Bonn.....	26
Bunge.....	24
Bard-Parker Handel.....	48
Bard-Parker Blade.....	48
Bishop-Harmon.....	50
Biploar Forceps.....	51
Biploar.....	51
Biplor Cord.....	51

## C

Colibri.....	01, 02
Castroviejo.....	12, 13, 14, 20, 30, 31, 36
Culler.....	18
Clayman.....	16, 19
Cook.....	35
Castroviejo.....	45
Chuck Handle.....	48
Cannulas.....	52, 53
Cattract Set.....	54

## D

Dyson-Colibri.....	01
Dressing.....	04
Desmarres.....	08
Drysdale.....	19
Drews.....	31, 49, 52

## E

Elsching.....	06, 20
Erhardt.....	11
Eye Scissors.....	17
Emery.....	42
Elliot.....	45
Enucleation Set.....	55, 59

## F

Farkas-Bracken.....	06
Fine.....	07

Francis.....	08, 43
Fixation Forceps.....	28
Fenzl.....	29
Feaster-Barraquer.....	34
Flieringe.....	39
Finskind-castroviejo.....	36

## G

Gill.....	01, 21, 43
Gills-Colibri.....	01
Girards.....	03
Gradle.....	07
Graefe.....	06
Gills-Vannas.....	15
Gimbel.....	18
Gills-Welsh.....	23
Girards.....	22
Gass.....	21
Graether.....	26, 24, 26, 47
Galand.....	28
Grandon-Barraquer.....	34
Guyton-Park.....	38
Golf-Club.....	43
Girard-swan.....	42

## H

Harms-Colibri.....	01
Hofmann-polack.....	02
Harms.....	07
Hunt.....	08
Hartman.....	10
Halsted.....	10
Hirschman.....	22, 27
Hydro-Dissection.....	49
Hydrodissection.....	50

## J

Jewelers.....	09
Jameson.....	11
Jaffe-Bechert.....	19
Jaffe.....	22, 30, 41
Jaffe-Knolle.....	26
Jaffe-Maltzman.....	28
Jaeger.....	41
J-Shaped Cannula.....	50
Jensen.....	49

## K

Kuhnt.....	06
Katena.....	07, 05
Katena-Vannas.....	15
Knolle.....	18, 26
Kraff.....	16
Kimura.....	20
Kelly.....	21
Kansas.....	23, 25
Keener-Arit.....	23
Keener.....	23
Kuglen.....	29, 38





# Index

Kalt.....31  
 Kartz-Barraquer.....35  
 Knapp.....40  
 Kelman.....44  
 Ketena.....45  
 Knolle-Kelman.....47  
 Kartz.....50,49  
 Kraff.....49  
 Knolle-Pearce.....52

## L

Lester.....04,23  
 Lambert.....08  
 Lordan.....08  
 Lindstrom.....13  
 Lewis.....20  
 Lieppman.....22,47  
 Lewcky.....27  
 Lancaster.....37  
 Lester-Burcs.....37  
 Lamellar Blade.....48  
 Lid set.....57  
 Lacrimal Set.....58

## M

Maumenee-Colibri.....03  
 Manhattan Eye & Ear.....05  
 McPherson-Vannas.....15  
 McPherson-Westcott.....14  
 Morrison.....23  
 Meyerhoefer.....21  
 McPherson.....30,36  
 Maltzman-Fenzl.....28  
 McInyre.....32  
 McIntyer.....34  
 Meaumennee-Park.....38  
 McIntrya.....47  
 Miniature Blade.....48  
 Mcintyre-Binkhorst.....52

## N

New Orleans.....20  
 Nightingale.....21  
 Nishi.....53

## O

Osher.....13

## P

Polack.....02  
 Pierse.....02  
 Pisacano.....18  
 Paton.....25

## R

Rubman Bores.....03  
 Rentsch"Boat Hook".....29  
 Reversed Simcoe.....46  
 Randolph.....44  
 Rubman-Katzin.....49

Radil Keratotomy Set.....56

## S

Sauer.....05  
 Suture Removal Forceps.....09  
 Serrefine.....10  
 Shepard II.....11  
 Shaepard.....16  
 Shepard.....26  
 Sinsky.....27  
 Sauer.....35  
 Smirmaul.....39  
 Schott.....39  
 Stevenson.....40  
 Stiletto Blade.....42  
 Spud & Gouge.....43  
 Simcoe.....46  
 Swiss.....46  
 Swiss Model.....45  
 Stiletto Knife.....48  
 Simcoe.....49  
 Silicone.....51  
 Sheet.....52  
 Simcoe.....52  
 Sterilizing Case.....58  
 Sterilizing Case.....59

## T

Tennant-Coilibri.....03  
 Tissue.....04  
 Tannant.....07,19  
 Tennant-Troutman.....06  
 Troutman.....06,33,44  
 Troutman-Castroviejo.....12  
 Twist-Hook.....28  
 Tennant.....35  
 Tooke.....43  
 Thurmond.....53

## V

Vannas.....15

## W

Wies.....06  
 Westcott.....14  
 Wilder.....20  
 Williams.....37  
 Wills Eye Hospital.....43  
 West.....44





Sheffield Surgical Building, Wazirabad Road, Ugoki, Sialkot-51310 Pakistan.  
Tell # +92-333-1400682  
Email: sheffieldsurgicalinstruments@gmail.com  
Email: info@sheffieldsurgicalinst.com  
Web: www.sheffieldsurgicalinst.com

